



ELLIPSOMETRY ELLIPSOMETER AND THIN FILMS

Dr. Gilles Benoit

Massachusetts Institute of Technology, Cambridge, MA

Dr. Richard Sun

Angstrom Sun Technologies Inc., Acton, MA



Outline



- ❖ Principle
- ❖ Instrumentation
- ❖ Modeling
- ❖ General Applications
- ❖ Summary



PRINCIPLE

On Spectroscopic Ellipsometry



Nature of the Light Waves



Light waves are electromagnetic in nature and require four basic field vectors to completely describe them.

- The electric-field strength: E
- The electric-displacement density: D
- The magnetic-field strength: H and
- The magnetic-flux density: B

Of these four vectors, the electric-field strength is chosen to describe the polarization state of light waves. It is because the force exerted on the electrons by E is much greater than by the magnetic-field of the wave when light interacts with matter. Once the polarization of electric-field vector has been determined, the other three vectors can be also determined based on Maxwell's field equations and the associated constitutive material relations.



Mathematic Description of Light Wave



For simplification, we consider a linearly polarized plane wave which propagates along the positive z direction and vibrates in the x-plane. The electric field of the wave can be written as:

$$c^2 \frac{\partial^2 E_x}{\partial z^2} = \varepsilon \frac{\partial^2 E_x}{\partial t^2} + 4\pi\sigma \frac{\partial E_x}{\partial t} \quad \Rightarrow \quad E_x = E_0 \exp\left[i\omega\left(t - \frac{z\tilde{N}}{c}\right)\right]$$

where E_x is the x-component of the electric-field strength; ε is the dielectric constant; σ is the conductivity.
 $\omega (=2\pi\nu)$ is the angular frequency, E_0 is the maximal value of the electric field strength.

\tilde{N} is the complex index of refraction

$$\tilde{N} = n - ik$$



$$\Rightarrow E_x = E_0 \exp\left[i\omega\left(t - \frac{z(n - ik)}{c}\right)\right] \Rightarrow E_x = E_0 \exp\left[-\frac{\omega k}{c} z\right] \cdot \exp\left[i\omega\left(t - \frac{zn}{c}\right)\right]$$

Damped Term

Undamped Term




Light Wave Propagation





Linearly and Circularly Polarized Light




The special cases of linear (a) and circular (b) polarization. In (a), the dashed line indicates the locus of the terminus of the electric vector E . In (b), L and R represent the left- and right-circular polarizations, respectively.



Elliptically Polarized Light





---- The extremity of its electric field vector describes an ellipse



(a) Elliptically polarized light; (b) decomposition of elliptically polarized light into two mutually perpendicular plane polarized waves P and S with a phase difference, Δ ;
R.W. Pohl, Optik und Atomphysik, Springer-Verlag, Berlin, 1958



Parameters to Describe the Ellipse




$$\sin 2\epsilon = \sin 2\psi \sin \Delta$$

$$\tan 2\theta = \tan 2\psi \cos \Delta$$




- (1) Azimuth θ , which is the angle between the major axis of the ellipse and the positive direction of the P. It defines the orientation of the ellipse ($-90^\circ \leq \theta < 90^\circ$).
- (2) Ellipticity e , which is the ratio of the length of the semi-minor axis and that of its semi-major axis. $e = a/b = \tan(\epsilon)$
- (3) Amplitude A , which is a measure of the strength of elliptical vibration. Its square is proportional to the energy density of the wave. $A = (a^2 + b^2)^{1/2}$
- (4) Handedness, which is used to describe wave propagation sense. Right handed – clockwise; Left handed – counter-clockwise when looking into the beam.

Plane of Incidence






P-polarization / S-polarization



Source: Newport



Poincaré Sphere



$$\sin 2\varepsilon = \sin 2\Psi \sin \Delta$$


$$\tan 2\theta = \tan 2\Psi \cos \Delta$$



Poincaré Sphere




R.M.A Azzam & N.M. Bashara: Ellipsometry and Polarized Light, Elsevier Sci., New York, 1999





Reflection of Light




$$\rho = \frac{R^P}{R^S} = \tan \psi \cdot e^{j\Delta} = f(n_i, k_i, d_i \dots)$$



Instrumentation




Typical SE Layout






Typical SE Layout





Spectroscopic Ellipsometers (SE)

Single Wavelength Ellipsometer (SWE)



VUV / PUV


Vis

Vis + NIR

IR-SE



Ellipsometer Classification





PRINCIPLE, INSTRUMENTATION AND MODELING

On Spectroscopic Ellipsometry

Modeling



Model and Its Analyses



$$\rho = \frac{R^P}{R^S} = \tan \psi \cdot e^{j\Delta} = f(n_i, k_i, d_i \dots)$$

Measured Data

Tan Ψ



Cos Δ



Physical Model

Estimated sample structure



- Film Stack and structure
- Material n, k, dispersion
- Composition Fraction of Mixture

Measurement
=
Calculation ?

Yes


T_i, n_i, k_i

No





Constants for BULK



2. Snell's Law:

$$\tilde{N}_a \sin \Phi_a = \tilde{N}_c \sin \Phi_c$$

3. Ellipsometry Equation

$$\rho = \frac{r_{ac}^P}{r_{ac}^S}$$

4. Optical Constants

$$r_{ac}^P = \frac{\tilde{N}_c \cos \Phi_a - \tilde{N}_a \cos \Phi_c}{\tilde{N}_c \cos \Phi_a + \tilde{N}_a \cos \Phi_c}$$

$$r_{ac}^S = \frac{\tilde{N}_a \cos \Phi_a - \tilde{N}_c \cos \Phi_c}{\tilde{N}_a \cos \Phi_a + \tilde{N}_c \cos \Phi_c}$$

$$\tilde{N} = \sin \Phi_a \sqrt{1 + \left(\frac{1 - \rho}{1 + \rho} \right)^2 \tan^2 \Phi_a}$$



Air/Film/Substrate




- Multi reflections between interface a/b and interface b/c
- From the Fresnel coefficients at the interface a/b and b/c, the total reflection coefficients can be calculated as:

$$R^P = \frac{r_{ab}^P + r_{bc}^P e^{-j2\beta}}{1 + r_{ab}^P r_{bc}^P e^{-j2\beta}}$$

$$R^S = \frac{r_{ab}^S + r_{bc}^S e^{-j2\beta}}{1 + r_{ab}^S r_{bc}^S e^{-j2\beta}}$$


$$\beta = 2\pi \left(\frac{d}{\lambda} \right) (n_b - jk_b) \cos \Phi_b$$



$$\rho = \frac{R^P}{R^S} = \tan \psi \cdot e^{j\Delta}$$



Analysis for Multilayer



- M films in structure
- Planar assumed
- Isotropic and homogenous
- Incident angle known

• Facts:

- M layers + one substrate
- M+1 interfaces

• Unknowns:

- Each layer thickness (m)
- Optical constants (n & k, 2m)

• Ellipsometer Data:

- One set of Ψ & Δ at each wavelength



Reflection and Transmission at Interface



1. Fresnel Reflection Coefficients:

$$r_{(i-1)/i}^P = \frac{\tilde{N}_i \cos \Phi_{i-1} - \tilde{N}_{i-1} \cos \Phi_i}{\tilde{N}_i \cos \Phi_{i-1} + \tilde{N}_{i-1} \cos \Phi_i}$$

$$r_{(i-1)/i}^S = \frac{\tilde{N}_{i-1} \cos \Phi_{i-1} - \tilde{N}_i \cos \Phi_i}{\tilde{N}_{i-1} \cos \Phi_{i-1} + \tilde{N}_i \cos \Phi_i}$$

2. Fresnel Transmission Coefficients:


$$t_{(i-1)/i}^P = \frac{2\tilde{N}_{i-1} \cos \Phi_{i-1}}{\tilde{N}_i \cos \Phi_{i-1} + \tilde{N}_{i-1} \cos \Phi_i}$$

$$t_{(i-1)/i}^S = \frac{2\tilde{N}_i \cos \Phi_i}{\tilde{N}_{i-1} \cos \Phi_{i-1} + \tilde{N}_i \cos \Phi_i}$$

3. Snell's Law $\tilde{N}_i \sin \Phi_i = \tilde{N}_{i+1} \sin \Phi_{i+1}$

4. Phase Factor

$$\beta_i = 2\pi \left(\frac{d_i}{\lambda} \right) \tilde{N}_i \cos \Phi_i$$





Ellipsometry Equations



Interface Matrix I (a & b)

$$I_{ab} = \begin{bmatrix} 1/t_{ab} & r_{ab}/t_{ab} \\ r_{ab}/t_{ab} & 1/t_{ab} \end{bmatrix} = (1/t_{ab}) \begin{bmatrix} 1 & r_{ab} \\ r_{ab} & 1 \end{bmatrix}$$

Layer Matrix


$$L = \begin{bmatrix} e^{j\beta} & 0 \\ 0 & e^{j\beta} \end{bmatrix}$$

$$S = I_{01} L_1 I_{12} L_2 I_{23} L_3 I_{34} L_4 \cdots I_{(m-1)m} L_m I_{ms} = \begin{bmatrix} S_{11} & S_{12} \\ S_{21} & S_{22} \end{bmatrix}$$

$$R = \frac{S_{21}}{S_{11}} \quad T = \frac{1}{S_{11}} \quad \longleftrightarrow \quad \rho = \frac{R^P}{R^S} = \tan \psi \cdot e^{j\Delta} = f(n_i, k_i, d_i \cdots)$$



Model and Its Analyses




$$\chi^2 = \frac{1}{2n - m - 1} \sum_{i=1}^n [(Tan\Psi_{Theory}^i - Tan\Psi_{Exp}^i)^2 + (Cos\Delta_{Theory}^i - Cos\Delta_{Exp}^i)^2]$$



How to Set up MODEL?



Layer Stack in SE Analysis






NK Table



	Wavelength(nm)	N	K
1	163.2	0.5251	2.1678
2	163.6	0.5329	2.1739
3	164.0	0.5403	2.1801
4	164.5	0.5473	2.1865
5	164.9	0.5538	2.1932
6	165.3	0.5600	2.2000
7	165.8	0.5658	2.2070
8	166.2	0.5713	2.2142
9	166.7	0.5765	2.2216
10	167.1	0.5814	2.2292
11	167.6	0.5861	2.2370
12	168.0	0.5905	2.2449
13	168.5	0.5947	2.2529
14	168.9	0.5987	2.2612
15	169.4	0.6025	2.2696
16	169.9	0.6062	2.2781
17	170.3	0.6098	2.2867
18	170.8	0.6133	2.2955
19	171.3	0.6167	2.3045
20	171.7	0.6200	2.3135
21	172.2	0.6234	2.3227
22	172.7	0.6266	2.3320
23	173.2	0.6300	2.3414
24	173.7	0.6333	2.3509
25	174.2	0.6367	2.3606
26	174.6	0.6403	2.3703
27	175.1	0.6439	2.3801
28	175.6	0.6477	2.3899
29	176.1	0.6516	2.3999
30	176.6	0.6559	2.4099
31	177.1	0.6605	2.4199
32	177.7	0.6656	2.4299
33	178.2	0.6711	2.4399
34	178.7	0.6769	2.4499


Optical Constants

◇ Si






Mixture of Materials






Mixture of Materials





Mixture of Materials (Interface)



a. Thickness
b. Volume Fraction of Lower Material

Properties
(Mixture of Above and Beneath Material)

NK Table:


- Homogeneous
- Graded

NK Profile:

- Linear
- Gaussian
- Exponential
-more


Physics Model:

- Bruggeman
- Maxwell-Garnett





Mixture of Materials (Surface Roughness)





Physics for Mixture of Materials

Microscopic Property

Polarizability: α

$$\frac{4\pi}{3} N\alpha = \frac{\epsilon - 1}{\epsilon + 2}$$

Macroscopic Property

Dielectric Response: ϵ

Clausius-Massotti Equation

(N is the density of polarizable units)



$$\frac{\epsilon - 1}{\epsilon + 2} = f_a \frac{\epsilon_a - 1}{\epsilon_a + 2} + f_b \frac{\epsilon_b - 1}{\epsilon_b + 2}$$

Lorentz-Lorentz Equation


(f's are volume fractions for each polarizable components a & b)



Maxwell-Garnett Equation

$$\frac{\epsilon - \epsilon_h}{\epsilon + 2\epsilon_h} = f_a \frac{\epsilon_a - \epsilon_h}{\epsilon_a + 2\epsilon_h} + f_b \frac{\epsilon_b - \epsilon_h}{\epsilon_b + 2\epsilon_h}$$

a and b are in a host with dielectric constant of ϵ_h
(Philos. Trans. R. Soc., London, 205, 1906, 237)







Effective Medium Approximation (EMA)

A.G. Bruggeman, Ann. Phys. (Liepzig) 24, 636 (1935)

The Bruggeman model is also called Effective Medium Approximation (**EMA**). In this model, the two media play exactly the same role. The effective dielectric function of the mixture is given by the second order equation :


$$\epsilon = \epsilon_h \quad \Rightarrow \quad 0 = f_a \frac{\epsilon_a - \epsilon}{\epsilon_a + 2\epsilon} + f_b \frac{\epsilon_b - \epsilon}{\epsilon_b + 2\epsilon}$$



The two materials play the same role but they are in interaction. Each type of inclusion is in interaction with the medium (it is supposed nevertheless spherical, depolarization coefficient equal to 1/3).




Surface and Interface






Direct and Non-Direct Technique






SE vs SEM for Inhomogeneous Film TiO_2



SE Model		SEM Photo
Thickness (nm)	<u>n@1550nm</u>	
13.5	1.7211	
424.42	2.2363	
66.71	2.2825	
Si Substrate		



Dispersion





Dispersion (UV Contribution)



Dispersion

UV Contribution

- Cauchy NK

$$n = A + B / \lambda^2 + C / \lambda^4$$

$$k = D + E / \lambda^2 + F / \lambda^4$$

- Cauchy $\varepsilon_r/\varepsilon_i$

$$\varepsilon_r = A + B / \lambda^2 + C / \lambda^4$$

$$\varepsilon_i = D / \lambda + E / \lambda^3 + F / \lambda^5$$

- Sellmeier

$$\varepsilon_r = 1 + (A^2 - 1) \lambda^2 / (\lambda^2 - B)$$

$$\varepsilon_i = C / \lambda + D / \lambda^2 + E / \lambda^3$$

- Exponential

$$\varepsilon_r = A + B \lambda^2 + C (D + (E + F / \lambda^2) / \lambda^2) / \lambda^2 \lambda^2$$

$$\varepsilon_i = 0$$

- Tauc-Lorentz

$$\varepsilon_i(E) = \frac{A E_0 C (E - E_g)^2}{((E^2 - E_0)^2 + C^2 E^2) E} \quad \text{for } E > E_g$$

$$\varepsilon_i(E) = 0$$

for $E < E_g$

- more



Dispersion (IR Contribution)

Dispersion



IR Drude

$$\varepsilon_r = P - (1/L_o)^2 \lambda^2 / (1 + ((1/\tau)\lambda)^2)$$

$$\varepsilon_i = (1/\tau)(1/L_o)^2 \lambda^3 / (1 + ((1/\tau)\lambda)^2)$$


Note: P is the polarization, $1/\tau$ the mean free pass and $1/L_o$ the inverse of the plasma wavelength.



Dispersion (Absorption Bands)



Dispersion



- Lorentz

$$\varepsilon_r = A\lambda^2(\lambda^2 - L_0^2) / [(\lambda^2 - L_0^2)^2 + \gamma^2\lambda^2]$$

$$\varepsilon_i = A\lambda^3\gamma / [(\lambda^2 - L_0^2)^2 + \gamma^2\lambda^2]$$

Note: A is the amplitude, L_0 is the central wavelength and γ is the width of the band or peak.

$$\varepsilon_r = A \cdot \text{imag}(\text{Double}W(z))$$

$$\varepsilon_i = A \cdot \text{real}(\text{Double}W(z))$$

$$\text{Double}W(z) = (\exp(-z^2) \cdot \text{erfc}(-iz))$$

$$z = (1/L_0 - 1/\lambda) / \gamma$$

- Gaussian



APPLICATIONS

with Spectroscopic Ellipsometry



Applications

- Optical constants or dielectric constants for bulk materials
- Optical constants and thickness for films
- Surface, Interface and Composites
- Alloy concentration determination
- Real time monitoring for growth or etching kinetics study
- Band gaps
- Porosity
- Other properties derived from optical /dielectric constants



General Applications



- **Materials:**


- Metals
- Polymers
- Ceramics and Glasses
- Semiconductors and its Compounds
- Composites

- **Application Fields**

- **Semiconductor Industry:** Photoresist, Gate dielectrics, Semiconductors and their alloys or compounds such as, SiGe, InGaAs, AlInGaAs
- **Photonics**
 - Optical coatings
 - Semiconductor compounds
 - Functional films in Optical MEMS
- **Data Storage**
 - Diamond-like carbon (DLC)
 - Magnetic films
- **Flat Panel Display (FPD)**
 - Thin film transistors (TFT) stack
 - Conductive oxide: Indium Tin Oxide (ITO)
- **Solar Cell Industry**




Applications







Multilayer Dielectric Stack



Tan(PSI)



Cos(Delta)





Summary



- Ellipsometry is a technique through monitoring the polarization state changes.
 - It offers nondestructive way to accurately determine film thickness and optical constants
- Many types of ellipsometers are available in the market. The selection and use of ellipsometer are highly application-oriented.
- SE has advantages over many other techniques.
- The ellipsometry is the MODEL based technique. It is not direct readout type instrument.
Therefore:
 - It is critical to use the right model to get right results



Ellipsometry Advantages



Ellipsometry has many advantages over other techniques:

1. More ***comprehensive*** than any other tools: Optical, electrical, physical (or structural) and chemical (composition, bonding). All these information may be obtained from only one measurement.
2. ***Nondestructive*** compared with SIMS for composition
3. ***Non-contact*** compared with four-point probe for conductivity
4. ***Non-contact*** and no pattern needed compared with stylus profilometer for thickness
5. ***No vacuum*** requirement compared with all e-beam or ion-beam based instrument



Advantages – cont.

6. Ellipsometry is an ***absolute*** technique (no need of reference or standards)
7. Ellipsometry gives ***twice*** more information (both phase and amplitude ratio) than reflectometry (only intensity). In addition, as ellipsometer measures the polarisation state and not the intensity, it is less sensitive to light intensity fluctuations.
8. The phase information from ellipsometry is very ***sensitive to surface*** layers. Therefore, it is the best non-destructive technique for thin film characterisation.



Further Information Available

- ▣ www.angstec.com
- ▣ Info@angstec.com
- ▣ Phone: (978) 263-6678
- ▣ Fax: (978) 263-6675