



Advanced Technique for Si_{1-x}Ge_x Characterization: Infrared Spectroscopic Ellipsometry

Richard Sun Angstrom Sun Technologies Inc., Acton, MA

Joint work with Darwin Enicks, I-Lih Teng, Janice Rubino ATMEL, Colorado Springs, CO

Presented at Advanced Semiconductor Manufacturing Conference 2004 Organized by SEMI

May 4-6, 2004, Boston



Outline



- 1. Introduction
- 2. Epilayer Growth
- 3. Methodology Characterization of SiGe films
- 4. Results and Discussion (IRSE vs. Xrd)
 - SiGe Box
 - Si Capped-box layers
 - Graded layers
 - Boron doped SiGe
 - Carbon alloyed SiGe with and without boron doped
- 5. Summary



Introduction



- Si_(1-x)Ge_x and Si_(1-x-y)Ge_xC_y films are top contenders for a wide range of high speed devices.
- Pseudomorphic growth of sub-50 nm SiGe films poses significant challenges for semiconductor manufacturers especially related to process development and control for enhanced manufacturability.
- Fast, reliable acquisition of in-line parametric data is critical to develop and sustain a manufacturable SiGe process.
 - %Ge, growth rate, dopant concentration, sheet resistance, strain effect
 - Across wafer uniformity
- Typical measurements are with Xrd, SIMS, and FPP which can be costly, "noisy", and time consuming for in-line qualifications.
- Spectroscopic ellipsometry is a non contact, non-destructive technique, which is suitable for both thin film and bulk characterization. Through accessing dielectric properties of the film, other properties can be extracted accordingly, such as physical, chemical and electrical properties.



Growth of Strained SiGe



- LPCVD (low pressure chemical vapor deposition) is performed in an AMAT 5200 Epi Centura Reactor.
- Pseudomorphic (lattice matched), sub-50 nm SiGe grown by thermal decomposition of silane (SiH₄), germane (GeH₄), methyl silane (CH₃SiH₃), and diborane (B₂H₆) in a H₂ ambient at ~600 °C.
 - SiH₄ + GeH₄ are Si and Ge sources for crystal growth
 - B_2H_6 is p-type dopant source for HBT base region.
 - CH₃SiH₃ is carbon source to minimize boron outdiffusion (captures Si interstitials) and provide strain compensation (reduces compressive strain)









Comparison between Different Techniques for SiGe **Technique**

Xrd

Strength

- repeatable, reliable, sub-20 nm capability,
- sensitive to lattice parameter variations,
- indicates strain relaxation. •
- non-destructive
- medium cost

Weakness

- time consuming,
- large spot size, not capable for product measurement.
- dopant incorporation may result in low Ge measurement due to lattice parameter reduction

FPP	 fast low cost good indicator of active dopant
SIMS	 sub-20 nm capability great indicator of dopant profiles and concentrations good measure of contaminants such as

oxygen

- difficulty for sub-50nm films
- influenced by dopants from adjacent layers
- not capable for product measurement, •
- destructive

- time consuming
- high cost
- does not indicate active or substitutional dopant incorporation
- destructive





Methodology - Characterization of SiGe

- UV-Visible Spectroscopic Ellipsometry (SE) measurement and analysis yield:
 - Thickness of SiGe layer and Si cap layer
 - Ge concentration in both box and graded SiGe layers Infrared SE (IR-SE) measurement and analysis yield:
 - Dopant concentration and electrical properties

 Both UV-Vis SE and IR-SE capabilities are integrated into IRSE tool (refer to presentation)

Model and Its Analyses

$$\rho = \frac{R^{P}}{R^{S}} = Tan\psi \cdot e^{j\Delta} = f(n_{i}, k_{i}, d_{i} \cdots)$$







Introduction to Alloy Model (2)

Transition Energy as a function of Ge content



Principle of Dopant Concentration Characterization (Drude Model)

$$\varepsilon_{1}(\omega) = \varepsilon_{\infty} - \frac{\omega_{p}^{2}}{\omega^{2} - \omega_{\tau}^{2}}$$
$$\varepsilon_{2}(\omega) = \frac{\omega_{p}^{2}\omega_{\tau}}{\omega(\omega^{2} - \omega_{\tau}^{2})}$$

 ω_p : plasma frequency

 ω_{τ} : scattering frequency

the conductivity of the material

$$\sigma = \varepsilon_0 \frac{\omega_p^2}{\omega_\tau}$$

the free carrier density N

$$\mathbf{N} = \boldsymbol{m}^* \frac{\boldsymbol{\varepsilon}_0 \boldsymbol{\omega}_p^2}{\mathbf{e}^2}$$

where m* is the effective mass

□ the free carrier mobility

$$\mu = \frac{\mathbf{e}}{\mathbf{m}^* \boldsymbol{\omega}_{\tau}}$$



Sample Description





Native Oxide ~180 Å Si Cap SiGe ~350Å Ge: ~25% + B or C Si Substrate **Native Oxide** ~180 Å Si Cap SiGeC ~350Å Ge: ~25% Si Substrate **Native Oxide** ~180 Å Si Cap SiGeC ~350Å Ge: ~25% + B

Si Substrate



Typical SE Spectra (box only)







Typical Spectra (with cap)











Sample	Ge Concentration				Thickness (Angstroms)						
	%Ge Box		%Ge Graded		SiGe Box		Graded SiGe		Si Cap		Nati∨e Oxide
ID.	IRSE	Norminal	IRSE	Norminal	IRSE	Norminal	IRSE	Norminal	IRSE	Norminal	IRSE
ww01017	22.2 ± 0.3	24.1	15.2	15.2	268.8	260	227.1	234			7.1 ± 0.2



Ge Concentration Profile



A MBIANT SIDE



Mapping on Si Cap Thickness



Presenter: Richard Sun (rsun@angstec.com)

Fab 5 Colorado Springs



Mapping on SiGe Layer Thickness

Fab 5 Colorado Springs





Mapping on Ge % in SiGe Box









Characterization – Summary Table

Sample	Ge Concentration				Thickness (Angstroms)						
	%Ge Box		%Ge Graded		SiGe Box		Graded SiGe		Si Cap		Nati∨e Oxide
ID	IRSE	Norminal	IRSE	Norminal	IRSE	Norminal	IRSE	Norminal	IRSE	Norminal	IRSE
ww01010	$\textbf{9.7} \pm \textbf{0.3}$	9.8			$\textbf{348.6} \pm \textbf{4.3}$	344					5.6 ± 0.2
ww01012	$\textbf{18.1} \pm \textbf{0.2}$	19.2			$\textbf{383.8} \pm \textbf{1.6}$	392					5.2 ± 0.2
ww01014	$\textbf{25.1} \pm \textbf{0.5}$	26.4			$\textbf{406.4} \pm \textbf{2.8}$	404					$\textbf{4.2} \pm \textbf{0.3}$
ww01009	9.7 ± 1.5	9.8			$\textbf{342.9} \pm \textbf{8.4}$	349			$\textbf{178.9} \pm \textbf{4.0}$	180	$\textbf{4.8} \pm \textbf{0.4}$
ww01011	19.0 ± 1.0	19.2			391.8 ± 13.9	388			$\textbf{178.9} \pm \textbf{3.2}$	181	$\textbf{4.9} \pm \textbf{0.6}$
ww01013	$\textbf{26.3} \pm \textbf{1.2}$	26.4			411.6 ± 12.6	407			$\textbf{178.5} \pm \textbf{2.6}$	175	4.6 ± 0.6
ww01015	$\textbf{22.2} \pm \textbf{0.3}$	24.2	15.3	15.6	256.3	251	237.6	236			7.7 ± 0.2
ww01016	$\textbf{21.9} \pm \textbf{0.3}$	24.0	15.2	15.3	273.3	254	225.3	235			3.9 ± 0.2
ww01017	$\textbf{22.2} \pm \textbf{0.3}$	24.1	15.2	15.2	268.8	260	227.1	234			7.1±0.2



IRSE - Spectra Comparison Doped vs. Non-doped Samples





WaveN.(cm-1)



Typical Fitting for Boron Doped Sample in Fab 5 Colorado Springs







Characterization Results from Infrared Channel

Sample	Resistance	Dopant Concentration		
(undoped 24-25 %)	na	na		
ww 3259 (plus boron)	228 ± 9 ohm/square	2.77 E+20 ± 4E+18		
ww 3260				
(plus carbon)	na	na		
ww 3261 (plus boron and carbon)	263 ± 15 ohm/square	2.71 E+20 ± 5 E+18		

Characterization Results from Visible Channel

Sample ID	SiGe Layer	Si Cap Thickness (Å)	SiGe Thickness (Å)	Ge %	IRSE - %Chg in Ge	Xrd - %Chg in Ge
Sample 1	SiGe only	171.7	389.5	24.9	0.00%	0.00%
Sample 2	SiGe + B	180.6	412.5	25.3	1.61%	-14.12%
Sample 3	SiGe + C	175.9	374.5	25.0	0.40%	-8.24%
Sample 4	SiGe + B + C	184.1	391.7	24.6	-1.20%	-18.61%



-0.8

2.0

2.5

3.0

Still Strained or Relaxed by Carbon Alloying?





3.5

4.0

4.5

Energy(eV)



- Strained model matches measured real SE Data.
- It is possible to tell whether SiGe layer is strained or relaxed with SE technique.







Precise measurements have been obtained for:

- Si cap layer thickness
- Si:Ge box layer thickness
- Ge concentration in Si:Ge box layer
- Ge gradient and thickness for graded Si:Ge layer
- Doped boron concentration and electrical properties







- IRSE method can reduce process qualification times significantly.
 - Typical Xrd "single point" measurements consume 1-2 hours/day of potential production time.
 - IRSE extracts single point measurement in minutes.
 - Typical Xrd "multipoint" measurements require 10 12 hours for 20+ points
 - IRSE extracts multipoint measurements in minutes
 - Weekly multipoint measurements are important for manufacturing control, but are difficult with Xrd due to time constraints.
 - Four Point Probe sheet resistance measurements on thin SiGe films are "noisy" and difficult due to probe contact problems, leakage to substrate, and probe conditioning issues.
 - IRSE is contact-less and offers improved reliability and speed.
- Fast, accurate determination of %Ge, thickness, sheet resistance, and dopant concentrations offers to increase production time, and enhance process development activity.





Further Information Available For any further information, please contact:

Email: info@angstec.com

www.angstec.com